





0

Mixing High Viscosity Solutions with the Vitl Co-Mix Laboratory Mixer

Abstract

The viscosity of a solution has a significant effect on the mixing efficiency of an orbital mixer. Viscous liquids such as glycerol and polyethylene-glycol (PEG) are widely used as a non-aqueous mixing solvent for drugs, vaccines and dyes, or as a stabilizing agent in protein solutions.

This work was carried out to examine the performance of the Vitl Co-Mix when used to mix solutions with a high concentration of glycerol in different microplates and tubes.

The results presented in this paper show that the Co-Mix is able to efficiently mix glycerol solutions of up to 90%wt (percentage by weight) in most vessel types.

Introduction

The Vitl Co-Mix combines a microplate orbital mixer and tube vortexer in a single unit, with independent speed controls for each module. It was developed to allow variable speed mixing of all types of multiple-well plates that are readily available and widely used in laboratories. The Co-Mix also comes with adapters for mixing an array of 0.5, 1.5 and 2ml microtubes.

Previous tests have shown that aqueous solutions of viscosity close to that of water are easily mixed within 1 minute (mostly less than 10 seconds) independent of the well geometry. However, it was decided to further examine the mixing capability of the Co-Mix when mixing highly viscous solutions, such as glycerol, which is widely used in the pharmaceutical and food industry.

The viscosity of a solution varies significantly with substance, concentration and temperature. The viscosity of glycerol reduces with increasing temperature and concentration. This is demonstrated in the References section at the end of this paper.

Materials and Method

The mixing performance of the Co-Mix was tested for different plate and tube types (vessels) with a range of fill volumes (from 10% to 100% of the maximum recommended working fill volume - WV) and at a tightly controlled ambient temperature of 20°C to 23°C.

The sample used was a solution of Glycerol (SIGMA-7893) in water at a concentration of 90%wt giving a viscosity of 219 cP at 20°C. The colourant was Ponceau 4R Red Dye at a ratio of 1 part in 30 by volume.

All tests were performed with the same viscosity solution - unless the solution failed to mix completely within one minute, in which case a lower viscosity solution (reduced concentration) was used and noted in the test results. The glycerol solution was dispensed into the sample vessel and the red dye pipetted to overlay it. Each combination was mixed on the Co-Mix for 60 seconds at the maximum permissible speed and the mixing efficiency was determined by visual inspection.



Application Note AN-001



Results and Conclusion

The tests were designed to examine the mixing performance of the Co-Mix with highly viscous solutions.

Table 1 below shows the tested vessel types, their physical properties and the achievable mixing viscosity for increasing fill volumes. As indicated by the test results, the vessel type and fill volume did not have a significant effect on the mixing efficiency other than for the 1.5ml micro-tubes.

The experiments demonstrated that the instrument is capable of mixing a 90%Wt glycerol solution, in most cases, independent of sample vessel geometry and fill volumes.

| Vessel Type | 96 Well MTP | 96 Well MTP | 96 Well PCR | 96 Well DWP | 384 Well MTP | 1.5ml Micro Tube | 2ml Micro Tube |
|---------------------------------|----------------|----------------|----------------|----------------|-----------------|------------------------------------------|----------------------|
| Vessel Shape | Circular | Circular | Circular | Circular | Square | Circular | Circular |
| Vessel Bottom | Flat | U-shaped | Conical | U-shaped | V-shaped | Conical | Round |
| Mixing Speed | 2000rpm | 2000rpm | 3000rpm | 3000rpm | 3000rpm | 2000rpm | 2000rpm |
| Vessel Volume | 350µl | 350µl | 150µl | 1200µl | 250µl | 1500µl | 2000µl |
| Working Volume (WV) | 100µl 28% | 100µl 28% | 100µl 66% | 500µl 42% | 100µl 40% | 1500µl 100% | 2000µl 100% |
| Fill Volume 10% of WV | 90%wt | 90%wt | 90%wt | 90%wt | 90%wt | 85%wt | 90%wt |
| Fill Volume 25% of WV | 90%wt | 90%wt | 90%wt | 90%wt | 90%wt | 85%wt | 90%wt |
| Fill Volume 50% of WV | 90%wt | 90%wt | 90%wt | 90%wt | 90%wt | 80%wt | 90%wt |
| Fill Volume 75% of WV | 90%wt | 90%wt | 90%wt | 90%wt | 90%wt | 80%wt | 80%wt |
| Fill Volume 100% of WV | 90%wt | 90%wt | 90%wt | 90%wt | 85%wt | Working Volume same as Maximum Volume | |
| | | | | | | | |
| Reference Viscosity | | 90%wt | | 85%wt | | 80%wt | |
| at 20°C | | 219cP | | 109cP | | 60cP | |

Table 1. vessel types, physical properties and achievable mixing viscosity for increasing fill volumes



Application Note AN-001



Pictures 1 and 2 below show the mixing of the 96 well flat-bottomed microtiter plate (MTP) at the five tested fill volumes of the glycerol and 1:30 dye solution:



Similarly, pictures 3 and 4 below show the mixing of the U-shaped 96 well deep well plate (DWP) at the five fill volumes:



It can be clearly seen from the images above that both sets of tests demonstrate excellent mixing within a short space of time without excessive vertical vortexing of the samples.



Application Note AN-001



Reference Data

The viscosity of glycerol at different concentrations (%wt percentage weight dilution in water) and temperatures is given in the table below. The viscosity data has been plotted for the three concentrations used in this experiment in the range 10°C to 50°C to allow the reader to interpolate from the quoted 20°C values.

Graph 1. Glycerol Viscosity /ºC

Table 2. Viscosity of Aqueous Glycerine Solutions

 in Centipoises/mPa s



Viscosity is specified in either the SI unit pascal-second (Pa.s) or in the CGS system (centimetregram-second) poise (p), or more commonly centipoise (cP). The relationship between these two

1 Pa.c = 1000 cP (or as per the table above) 1cp = 1mPa.s



units of measurement is: