


Behind The Scenes

– Verification Method: Dropcount




Behind the scenes

VERIFICATION METHOD: DROPCOUNT

After creating a liquid class or after assigning a liquid class to your sample liquid, the liquid class needs to be verified. Therefore, our verification method 'the Dropcount' was developed.

The idea behind it is that you fill the well with a defined volume (here 22µl) and you are dispensing spots with one per cent of the filling level (here 220 nL). If the liquid class fits to the respective liquid, there would be exactly 100 spots dispensed – if not, the offset in percent is directly correlating to the number of missing or additional spots.



1.

Assign the chosen liquid class to the liquid named "Dropcount"


Open Assay Studio, navigate to the liquid library and edit the liquid named "Dropcount". Assign the chosen liquid class to this liquid.

2.

Open the Dropcount protocol

Open the Dropcount protocol which is saved on your device. The protocol shows 30 columns of 4 target spots.

The goal of the verification run is to get 25 columns ($25 \times 4 = 100$) dispensed successfully. The selection of 30 columns of target spots just allows you to detect an underdispensing offset of the liquid class (more than 100 spots dispensed).



3.

Fill the source well (here with 22 µl)

Take your sample liquid and fill your source well. In this step, it is crucial that the pipette uses is calibrated and no pipetting error is introduced.

4.


Run the protocol

Click on "Dispense" and wait until the protocol is dispensed

5.

Read out the successfully dispensed spots

– compare the results in the software and on the plate.



Interpretation of the results

Too few droplets:

- One spot is containing more than 220 nL
- The droplets at the used dosing energy are bigger as deposited in the liquid class
- That's why the well is running empty too early
- The device is OVERDISPENSING

Action points:

Either:

Increase the droplet volume in the liquid class according to the offset in percent

Or:

Try the next lower viscosity liquid class

Too many droplets:

- One spot is smaller than 220 nL
- The droplets at the used dosing energy are bigger as deposited in the liquid class
- That's why the well is running empty too late
- The device is UNDERDISPENSING

Action points:

Either:

Decrease the droplet volume in the liquid class according to the offset in percent

Or:

Try the next higher viscosity liquid class



©2021 BICO AB. All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of BICO is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. BICO provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. BICO may refer to the products or services offered by other companies by their brand name or company name solely for clarity and does not claim any rights to those third-party marks or names. BICO products may be covered by one or more patents. The use of products described herein is subject to BICO's terms and conditions of sale and such other terms that have been agreed to in writing between BICO and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of BICO products in practicing the methods set forth herein has not been validated by BICO, and such nonvalidated use is NOT COVERED BY BICO'S STANDARD WARRANTY, AND BICO HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner BICO's terms and conditions of sale for the instruments, consumables or software mentioned, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by BICO that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact

Tel: +49 (0) 711 490 544 00

Email: info@dispendix.com

Website: www.dispendix.com