


Dispensing Detergents With the I.DOT


Elisa Kullick, Dispendix GmbH



Dispensing Detergents With The I.DOT

Detergents, also termed as tensides or surfactants are used in laboratories every day. They are commonly used in biochemistry, cell biology or molecular biology. A few examples of applications include cell lysis, protein solubilization, protein crystallization, or reduction of background staining in blotting experiments.

Detergents are amphipathic compounds with both hydrophobic and hydrophilic sides within the same molecule.



Liquid Class Identification & Dispensing Of Detergents

The liquid handling of detergents is challenging, as they tend to foam easily. Additionally, due to their ability to decrease the surface tension of water, detergent dilutions up to a specific percentage are often more fluid than pure water.

However, they start to get more viscous than water when the amount of detergent reaches an amount of more than ~ 5-10 %.

With the right settings and the identification of the fitting liquid classes the I.DOT makes the dispensing of detergents very simple. We could successfully dispense the following detergents with the I.DOT:

Liquid	Liquid Class	Dosing Energy Range	Droplet Volume Range	Well Type
Brij35 0,6 %	Water LQC	70 – 200 mbar*ms	13,5-33,1 nL	S.100
DDM (N-Dodecyl-B-D-Maltoside) 0,6 %	LQC0.8	70 – 200 mbar*ms	15,4-37,7 nL	S.100
Triton X-100 0,1 %	Water LQC	70 – 200 mbar*ms	13,5-33,1 nL	S.100
Triton X-100 1 %	LQC 0.7	70 – 200 mbar*ms	17,5-40,5 nL	S.100
Triton X-100 10 %	LQC2	70 – 200 mbar*ms	10,8-27,0 nL	S.100
TWEEN 20 0,1 %	Water LQC	70 – 200 mbar*ms	13,5-33,1 nL	S.100
TWEEN 20 1 %	Water LQC	70 – 200 mbar*ms	13,5-33,1 nL	S.100
TWEEN 20 10 %	LQC1	70 – 200 mbar*ms	12,3-28,8 nL	S.100
IGEPAL 0,1 %	Water LQC	70 – 200 mbar*ms	13,5-33,1 nL	S.100
IGEPAL 1 %	Water LQC	70 – 200 mbar*ms	13,5-33,1 nL	S.100
IGEPAL 10 %	LQC1	70 – 200 mbar*ms	12,3-28,8 nL	S.100



©2021 BICO AB. All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of BICO is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. BICO provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. BICO may refer to the products or services offered by other companies by their brand name or company name solely for clarity and does not claim any rights to those third-party marks or names. BICO products may be covered by one or more patents. The use of products described herein is subject to BICO's terms and conditions of sale and such other terms that have been agreed to in writing between BICO and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The use of BICO products in practicing the methods set forth herein has not been validated by BICO, and such nonvalidated use is NOT COVERED BY BICO'S STANDARD WARRANTY, AND BICO HEREBY DISCLAIMS ANY AND ALL WARRANTIES FOR SUCH USE. Nothing in this document should be construed as altering, waiving or amending in any manner BICO's terms and conditions of sale for the instruments, consumables or software mentioned, including without limitation such terms and conditions relating to certain use restrictions, limited license, warranty and limitation of liability, and nothing in this document shall be deemed to be Documentation, as that term is set forth in such terms and conditions of sale. Nothing in this document shall be construed as any representation by BICO that it currently or will at any time in the future offer or in any way support any application set forth herein.

Contact

Tel: +49 (0) 711 490 544 00

Email: info@dispendix.comWebsite: www.dispendix.com