

CAPP

D E N M A R K

CappTronic

ELECTRONIC PIPETTES

capptronic.capp.dk

THE CAPPTRONIC™ IS CERTAINLY ONE OF THE EASIEST ELECTRONIC PIPETTES TO SET AND USE. HAVING ONLY A ONE-SIDE PROGRAMMING MANUAL, FAMILIARIZING ONE'S SELF WITH CAPPTRONIC'S™ FUNCTIONS IS A MATTER OF MINUTES. APART FROM STANDARD PIPETTING, YOU CAN CHOOSE BETWEEN THE MOST WIDELY USED PIPETTING MODES INCLUDING REVERSE PIPETTING, MULTI DISPENSING, MIXING AND SEQUENTIAL DISPENSING. UP TO NINE INDIVIDUAL PROTOCOLS CAN BE STORED THANKS TO THE EXTENSIVE MEMORY CAPACITY. SUPPLIED WITH TWO LONG-LIFE AND EASILY EXCHANGEABLE LITHIUM BATTERIES, THE CAPPTRONIC™ ALLOWS FOR UNINTERRUPTED PIPETTING WITHOUT THE NEED TO STOP WORK IN ORDER TO RECHARGE.

THE O-RINGS ON THE TIP BASES PROVIDE AN EXCELLENT FIT WITH A GREAT VARIETY OF UNIVERSAL FITTING TIPS.



 CAPP APS
LANDBRUGSVEJ 10
DK-5260 ODENSE
DENMARK

TEL. (+45) 6613 6140
FAX. (+45) 6613 2770

MAIL: INFO@CAPP.DK
WEB: WWW.CAPP.DK



Ordering Info

Capp Tronic Single

CAT. No.	DESCRIPTION	INACCURACY %	IMPRECISION %	COLOR
E20-1BA	CAPP ELECTRONIC PIPETTE, SINGLE CHANNEL, 2-20 UL	1.00	0.30	
E200-1BA	CAPP ELECTRONIC PIPETTE, SINGLE CHANNEL, 10-200 UL	0.60	0.15	
E1000-1BA	CAPP ELECTRONIC PIPETTE, SINGLE CHANNEL, 100-1000 UL	0.50	0.15	

Capp Tronic Multi

CAT. No.	DESCRIPTION	INACCURACY %	IMPRECISION %	COLOR
E20-8BA	CAPP ELECTRONIC PIPETTE, 8-CHANNEL, 2-20 UL	1.00	0.30	
E200-8BA	CAPP ELECTRONIC PIPETTE, 8-CHANNEL, 10-200 UL	0.60	0.15	
E1200-8BA	CAPP ELECTRONIC PIPETTE, 8-CHANNEL, 100-1200 UL	0.50	0.15	

* ALL ELECTRONIC PIPETTES INCLUDE 2 LITHIUM BATTERIES