

we deliver  
you discover



THREE  
LEVELS OF  
SERVICE  
AVAILABLE!

## SERVICE PLANS FOR CS700-SERIES SCREW CAP RECAPPER

Four plans available to fit your needs.

1A



### STANDARD PM SERVICE PLAN

*Ship unit to workshop*

- One (1) scheduled PM inspection (recommended every 12 months)
- Performance test
- Thorough cleaning and lubrication of machine axis, spindles and guidings
- Functionality testing of machine axis, motor, push-off unit, optional tube detection sensor, slip couplings, bits and o-rings
- Inspection of screw head, fuses, guides, o-rings, springs, power cords and plugs (adjustment, replacement or repair of these parts done as necessary)
- Exclusions: all repair parts

1B




### STANDARD PM SERVICE PLAN

*Repair unit on-site*

- One (1) scheduled PM inspection (recommended every 12 months)
- Travel expenses related to PM visit
- Performance test
- Thorough cleaning and lubrication of machine axis, spindles and guidings
- Functionality testing of machine axis, motor, push-off unit, optional tube detection sensor, slip couplings, bits and o-rings
- Inspection of screw head, fuses, guides, o-rings, springs, power cords and plugs (adjustment, replacement or repair of these parts done as necessary)
- Exclusions: all repair parts

2




### LEVEL 2 SERVICE CONTRACT

*Repair unit on-site*

- One (1) annual PM inspection (includes all offerings in Standard PM)
- One emergency, on-site service visit
- Service visit includes repair parts\*, labor cost, and travel
- \*Exclusions: screw head

3



### LEVEL 3 SERVICE CONTRACT

*Repair unit on-site*

- Includes all offerings in Level Two
- Extra emergency on-site service visit
- Micronic software upgrades
- Service visit includes repair parts\*, labor cost and travel
- \*Exclusions: screw head

Need a quote? Have a question? Contact us today:



+31 (0)36 549 1010



info@nbsscientific.nl



www.nbsscientific.nl